


## Cardiovascular Guide Catheter (PTFE lined)

Precision, Performance, Perfection.



- **Functional lumen sizes:** 5FR .056", 6FR .070", 7FR .080", 8FR .090"
- **Variable shaft stiffness:** Ease of placement, exceptional pushability, torque control and product support
- **High tensile braid:** Improved kink resistance, increased pushability, improved torque
- **PTFE Liner:** exceptional inner surface lubricity
- **Manufacturing process:** providing reliability, super smooth surface and ease of use


### Guide Catheter Specification

FR Size	5FR	6FR	7FR	8FR
ID +/-0.001"	0.056"	0.0700"	0.0800"	0.0900"
OD +/-0.002"	0.0680"	0.0790"	0.0930"	0.1040"
Cushion Tip Length +1/-0mm	2 mm	2 mm	2 mm	2 mm
PTFE Lined Braided Shaft +/-2.5cm	3 sections	3 sections	2 sections	2 sections
Assembled Shaft Length +/-2cm	104 cm	104 cm	104 cm	104 cm

### Catheter Curve Setting Technology:

Comes from years of experience and multiple doctors' feedback

- \* 5 and 6F have 4 sections, 7 and 8F have 3 sections
- \* Proprietary Processing
- \* Detailed drawings and curve forms for 48 curves for 5, 6, 7, and 8F
- \* Perfusion ports as required



Contact us:

12 Glens Falls Technical Park, Glens Falls, New York 12801 USA

**E-mail:** [sales@precisionextrusion.com](mailto:sales@precisionextrusion.com) **Web Site:** [www.precisionextrusion.com](http://www.precisionextrusion.com) **Phone:** +1 (877)-792-1199