

CATHETER Solutions

- Angiographic
- Guiding
- Micro Catheter
- Electrophysiology



Capabilities

- Braiding
- Welding
- Reflowing
- RF Tip Forming
- Hole Punching
- Tapered

PVC PET COMPOUND PEEK PA PUR PC PE PP

Your Partner In EXCELLENCE

Precision Extrusion, Inc. of Glens Falls, New York, is a leader in the design and manufacturing of high quality thermoplastic tubing for the medical industry.

We are committed to unsurpassed quality control, ongoing investments in state-of-the-art manufacturing equipment and customer service that is second-to-none. Our highly skilled, talented team is dedicated to being your partner in excellence.

At Precision Extrusion, Inc., precision, performance and perfection really do matter.

PRECISION PERFORMANCE PERFECTION



All Content Copyright © 1995 - 2013 Precision Extrusion Inc. All rights reserved.



Precision: Numbers Speak For Themselves

Solutions by the numbers. Precision Extrusion Inc. is a company that measures success in small steps. We know how important each component of a project is, especially for small catheters. Our engineering team can assist in the development of custom tubing solutions, in any number of resins, to help you develop custom catheter.

- Most engineering plastics, compounds, and color matches
- Minimum diameter: 0.0025"; minimum wall thickness: 0.002"
- Minimum tolerance: 0.0005"

Performance: Precise Products, Exceptional Service

Solutions by the process. We put problems under a microscope. We custom engineer tubing solutions to the most demanding standards and work closely with our clients to develop a product that performs exactly as it should.

- ISO Certified
- Precise computerized equipment, constant screw RPM, inline lasers monitor and control OD, SPC for ID and wall.
- Prototyping and in-house tooling design
- Short lead times down to one week

Perfection: Single Lumen, Multi Lumen, Braided Catheter, Acetal Mandrel

Solutions by Engineering. At PEI, our skilled engineers have the ability to create unique solutions to tubing problems with the closest of tolerances. Let our team help create your next great device - on time, on budget, and to the most exacting standards. The word precision is right in our name. Let us offer our expertise on your next project.

- Engineering design consultation
- Complex configuration - coextruded, bump, twin, coiled, striped, over jacketed, and radiopaque tubing
- Secondary operations - tipping, flaring, RF welding, grinding, reflowing, precision cutting, and more.



12 Glens Falls Technical Park
Glens Falls, NY 12001
Phone +1(518) 792-1199
Toll Free +1(877) 792-1199
Fax +1(518) 792-3805

www.precisionextrusion.com

ISO Certified

"Engineers appreciate good engineering"

